

# Guided Surgery Workflow



## THE CUSTOMER SENDS TO THE MILLING CENTER:

- CT scan in DICOM format
- Duplicate diagnostic wax-up
- Models
- Medical Prescription
- Evobite



## THE MILLING CENTER PREPARES:

- Scanning Models with Evobite
- Design with software
- Fabrication of the surgical template
- Possible immediate-load Provisional



## THE CUSTOMER RECEIVES THE FINAL WORK:

- CT scan in DICOM format
- Duplicate diagnostic wax-up
- Models
- Medical Prescription
- Evobite
- +
- Surgical template
- Chewing index
- 3D printed model
- Design summary
- Eventual immediately loaded temporary